

The resuscitator with improved recoil time

The Ambu SPUR II provides users with exceptional tactile and visual feedback during resuscitation. The bag is highly responsive, with minimal mechanical resistance – some say it's a natural extension of the hand. The bag itself is thin and responsive and our SafeGrip™ surface ensures that the Ambu SPUR II can be firmly squeezed.

Our convenient handle makes it easier hold the unit and to ventilate a patient with just one hand if necessary. And our characteristic design provides optimum stroke volume with perfect recoil – yet is smaller and lighter than most other disposable units.

Optimal control

Our unique single-shutter valve has proven its reliability over two decades of emergency and

hospital use. The crystal clear valve housing provides an unimpeded view of the valve operation, and the slim design makes it easier to view the mask.

Adult and pediatric sizes have a soft splashguard for user safety and all pediatric and infant sizes are equipped with an override valve system that helps prevent too high lung pressure. All bags are colour-coded for fast identification of infant, pediatric and adult sizes.

The best resuscitator you'll ever throw away

Ambu SPUR II is the only resuscitator that is made from a SEBS material instead of PVC. This classifies Ambu SPUR II as environmentally safe and disposable, thus eliminating all risks of cross contamination. Simple is good. But simpler is better.

Key features of Ambu® SPUR® II

- Unique single-shutter valve system for reliable functionality
- Integrated handle for user comfort and uniform compression
- Thin-walled compression bag allows for lung compliance and "feel"
- SafeGrip™ surface for secure handling in stressful environments
- Swivel between valve and mask permits 360° positioning of the Ambu® SPUR® II in relation to the patient
- Easy attachment of manometer and PEEP valve



Ambu® SPUR® II

The unmistakable feeling of superior control

US: RX only

Ambu A/S
Baltorbakken 13
DK - 2750 Ballerup
Tlf.: +45 7225 2000
Mail: hello@ambu.com
www.ambu.com

Ambu® SPUR® II

The unmistakable feeling of superior control

- Improved recoil time and superior control
- No risk of cross contamination
- Environmentally safe and disposable



Our work

Building on heritage to invent the future

Bringing our finest moulding experts together, we created Spur II. A resuscitator with an optimised shape in a PVC-free material, that ensures the world's best recoil time and superior control. Ease of handling and no risk of cross contamination is what makes the difference at the end of the day.



Your Idea

Can you recreate "the Ambu feel" in a disposable device?

From hands-on experience, you are most likely familiar with the Ambu feeling: The unmistakable feeling of precision that has characterised our high quality resuscitator, since we invented the world's first self-inflating resuscitator more than 50 years ago. When asked if you felt there was room for improvement, your call was clear: Could we add the Ambu feeling to a disposable resuscitator.



Everyone's life

Designed to be trusted. Built for action

8 million uses over the past five years underline that a disposable and ready-to-use resuscitator deserves a place in every ambulance and across the hospital. Resuscitation should never be delayed and always be as precise as possible. Every second counts.

