

King Vision®

Video laryngoscope



Key Benefits

Portable

- Disposable blades allow economical use of the King Vision® video laryngoscope for all intubations.
- With the King Vision® you can improve patient care by having high performance visualisation capabilities for every intubation, at an affordable price.

Affordable

- The King Vision® is light weight, self-contained, battery operated and water resistant.
- The King Vision® can be used continuously about 80 minutes making it highly portable.
- Blades are individually packaged. Each blade is disposable, eliminating concerns on cross contamination.

Durable

- Designed to be your primary tool for all intubations.
- Video output capability is compatible with external monitor and recording devices.
- The robust, full-colour, non-glare display can resist repeated cleaning and normal use wear and tear.

King Vision® video laryngoscope

Visionary performance in the palm of your hand. We know that reliable performance is a necessity, not an option. When it comes to visualising the airway for routine and difficult intubations, you need a tool that you can depend on each and every time. The King Vision® video laryngoscope is designed with your expectations in mind.

Portable


The King Vision® display is light weight, self-contained, battery operated, and water resistant. The reusable display comes packaged in a protective, foam case. Blades are individually packaged making the King Vision® very portable.

Affordable

The disposable blades allow economical use of the King Vision® for all of your intubations. The low cost per procedure combined with the high performance visualisation capabilities make this the perfect value for the best patient care. Take the guess-work out of your intubations. See the vocal cords for every placement. Use the King Vision® for all of your intubations.

Durable

The King Vision® is designed to be your primary tool for intubations. The high-quality, made in America, display comes with a 1-year warranty. The robust, full-colour, smudge-proof, non-glare display can resist repeated cleaning and normal use wear and tear.



Ordering specifications

Item no.	Description	Quantity
KVL03	Standard Blade, Size 3	10
KVL03C	Channeled Blade, Size 3	10
KVIS01	Reusable Display	1
KVLKIT3	Vision Display, Case, AAA Batteries, 3 Channeled and 1 Standard Size 3 Vision Blades	1
KVCABL	Custom Video-out Cable (2.75 m / 9')	1

Sales Company
Ambu A/S
 Baltorpbakken 13
 DK - 2750 Ballerup
 Denmark
 Tel. +45 7225 2000
 Fax +45 7225 2053
 www.ambu.com

Manufacturer
King Systems
 15011 Herriman Boulevard
 Noblesville, IN 46060
 Tel. 317 776 6823
 Tel. 800 642 5464
 Fax 317 776 6827
 www.Kingsystems.com

CE US: Rx only