

Ambu® ResCue Mask



Key Benefits

- The ResCue Mask is transparent to allow visual check of patient's condition
- Head strap for fixation and secure positioning.
- For use on both adults and children
- Reusable facemask that is easy to clean, which extends the product life
- Available in softcase
- Phthalate-free material

Ambu® ResCue Mask

The Ambu ResCue Mask is designed to provide mouth-to-mask ventilation.

As a physical barrier between the rescuer and the patient Ambu ResCue Mask is used to give protected rescue breathing for any patient (any age group) and may also be used with a resuscitator.

The Ambu ResCue Mask helps protect rescuers in an emergency situation. It prevents direct contact with patient's mouth, nose and face and helps to overcome hesitation to start resuscitation.

Specifications

Product	
Size	
Dimensions	138 mm x 113 mm x 50 mm
Weight	Approximately 110g
Materials	
Mask	Polyvinylchloride
Softcase	Polyurethane, Ethylene Vinyl Acetate and Micro Fabric

Ordering specifications

Item no.	Ambu Rescue Mask
252 301 004	Ambu ResCue Mask – Adult mask with oxygen inlet, with elastic head strap and one-way valve with bacterial filter. Supplied in red softcase
252 301 104	Ambu ResCue Mask – Adult mask without oxygen inlet, with elastic head strap and one-way valve with bacterial filter. Supplied in red softcase

*ResCue Mask with Oxygen Inlet is Rx Only in the USA