

Six good reasons to use aScope™ 3

- aScope™ 3 offers improved patient safety as it is always available when needed
- aScope™ 3 is sterile with no risk of cross-contamination
- aScope™ 3 is portable, intuitive and fast to set-up and use
- aScope™ 3 gives you 3 sizes in one system at no additional cost
- aScope™ 3 is cost effective. Low investment costs. No repair and reprocessing costs
- aScope™ 3 gives you easy access to practice your bronchoscopy skills



Ambu® aView™
Portable, intuitive monitor for
Ambu® aScope™ 3 Family

Ambu® aView™ is a portable high resolution monitor that is compatible with all of Ambu's single-use flexible scopes in the Ambu® aScope™ 3 Family.

The clear and high resolution image quality enables easy navigation and fast identification of anatomical landmarks.

Key benefits of aView™ include

- Portable touch-screen monitor
- Operational in seconds
- Intuitive user interface
- Easy file transfer by USB
- Minimum 3 hours battery time
- Enhanced viewing with Image Extend
- Display of live image on external monitor, e.g. for training



Image Extend enlarges the live image on aView™ for better viewing at a distance.



Video-out functionality allows you to display the live image of aView™ on an external monitor, e.g. for training.

		SVF Placement of DLTs and BBs	Difficult Intubation	Bronchial Inspection	BAL/BW	Protected specimen brush sampling	PDT	Management of retained secretion	Management of tough retained secretion	
	Ambu® aScope™ 3 Slim 3.8/1.2 Ref # 402001000	●	●	●						ETT ≥ 5 DLTs ≥ 37
	Ambu® aScope™ 3 Regular 5.0/2.2 Ref # 403001000	●	●	●	●	●	●	●		ETT ≥ 6 DLTs ≥ 39
	Ambu® aScope™ 3 Large 5.8/2.8 Ref # 414001000		●	●	●	●	●	●	●	ETT ≥ 7
	Ambu® aView™ Ref # 405002000									

The aScope™ 3 Family consists of 3 bronchoscope sizes. The colour coding makes it easy and fast to find the right size.

See for yourself at visualisation.ambu.com

Ambu A/S
Baltorpbakken 13
DK-2750 Ballerup
Denmark
Tel.: 1 410 768 6464
Fax: 1 800 262 8673
www.ambuUSA.com



496 3605 01 Vol 1 2016/03



Ambu® aScope™ 3 Family
Tried and trusted

Ambu® aScope™ 3 Family ensures
flexible bronchoscopes are always
readily available



Better patient safety

aScope™ 3 ensures immediate access to a sterile flexible scope and eliminates the risk of cross-contamination, so you and your patients can feel safe.

Improved productivity

The single-use aScope™ 3 helps you save time and work smarter by eliminating time consuming steps required today to use, maintain and handle a reusable scope.

One system with three bronchoscope sizes to cover your needs

Whether you need a bronchoscope for intubation, placement of a DLT or suctioning of heavy secretion, the aScope™ 3 Family covers your needs. The colour coding makes it easy and fast to find the right size.



With the aScope™ 3 you always have a sterile bronchoscope available whenever you need it

aScope™ 3 is the safe and cost effective choice for fibreoptic intubation in the OR and ER



Improving patient safety

You and your patients can feel safe as the aScope™ 3 ensures you always have a flexible scope readily available whenever needed.

Procedures made simple

The sterile, single-use and ready-to-go concept of aScope™ 3 simplifies flexible intubation. The aScope™ 3 system is highly portable, intuitive and fast to set up. After use there is no handling or reprocessing; simply discard the scope.


Cost effective

aScope™ 3 offers low investment costs and no repair and reprocessing costs. And you can stop worrying about surgery delays due to lack of an available flexible scope.



aScope™ 3 is ideal for a wide range of procedures in the OR:

- Intubation
- Airway inspection
- Placing and checking of Double Lumen Tubes
- Placement of Bronchial Blockers



aScope™ 3 makes lack of available bronchoscopes a thing of the past in the ICU



Designed for bronchoscopy in ICU

Whether you need to intubate, perform a PDT or BAL procedure, or require high suction capacity for tough retained secretion, the three scope sizes in the aScope™ 3 Family cover your needs. The unique colour coding of the scopes makes it easy to select the right size. aScope™ 3 is immediate available, portable, and easy to use, making it the ideal choice for bedside bronchoscopy.

Your patient is safe

With the sterile single-use concept of aScope™ 3 your patient is safe as there is no risk of cross-contamination. And the ready-to-go concept ensures that you always have a scope readily at hand when you need it.

Cost effective

You can treat your patient here and now without having to wait for a scope to be found and brought to you. And you no longer have to worry about your scopes being broken and expensive repair costs.

aScope™ 3 is ideal for a wide range of bedside procedures in the ICU:

- Airway inspection
- Difficult intubation
- Management of retained secretion
- Bronchial Aveolar Lavage
- Percutaneous Dilatational Tracheostomy
- Protected Specimen Brush sampling

